

Clinical Cases

Localized area : Canine with dentin hyper sensitivity



Clean the surface



Apply



After 30 sec. application and rinsed

Large area : After cleaning



Polish after scaling
(mechanical instrumentation)



Dispense Gel Desensitizer
to a rubber cup




Rub each tooth surface with
low to moderate speed for 5 sec.

Pictures of product



Contents of the kit

- ① Gel Desensitizer
1 (3ml)
- ② 10 Plastic Needle Tip
+ 1 grey Needle Cap



Sold Separately

- 50 Plastic Needle Tip
+ 2 grey Needle Caps

Made in Japan by

SUN MEDICAL CO.,LTD.

571-2 Furutaka-cho, Moriyama,
Shiga, 524-0044, Japan
Phone:81-77-582-9981 Fax:81-77-582-9984
<http://www.sunmedical.co.jp>

Read carefully instructions and labels before the use - Notify body: TÜV SÜD Product Service GmbH - APRIL 2019 version

Exclusive European distributor :


GENERIQUE INTERNATIONAL
77/79 rue Pierre Brossolette • 93160 NOISY LE GRAND • Tel. 01 43 03 06 84
E-mail : generinter@aol.com - Web : www.generiqueinternational.com



Gel Desensitizer

Desensitizer for Hypersensitive Teeth

New gel-type desensitizer
Easier application and higher effectiveness



Easier application and higher effectiveness! New gel-type desensitizer

Easy application without
the hassles of multiple steps!
No measuring, mixing,
applicators or light-curing.



Point 1 Direct application from the syringe

Point 2 No need of rubbing motion

Point 3 Perfect viscosity lets it stay on the spot

Point 4 Biologically accepted for gingiva

Clinical procedure

1. CLEAN

Clean a tooth surface with
brushes or cotton pledgets.

2. APPLY

Form a thin layer and wait for 30 seconds.

The gel has a pleasant sweet flavor.

3.

**ASK
THE PATIENT
TO RINSE
WITH WATER**

TIP 1

If the area is too sensi-
sitive to have deposits
removed, you can ap-
ply the gel over the de-
posits with cotton ap-
plicator and rub it in.

TIP 2

Use a rotary cup with low
speed for a larger area.


Excellent performance on desensitization

MS Polymer

Oxalic acid

Potassium salt

Gel Desensitizer inherited MS Coat series' excellent performance in managing hypersensitivity, advanced with additional potassium salt. Nano sized MS polymer and oxalic acid react chemically with calcium in teeth and form a protective film containing fluoride (sodium fluoride) and potassium salt. The gel remains longer on teeth so dentin tubules can be sealed tighter than other MS Coat series




Increased acid resistance

Sodium fluoride 900 ppm


Included sodium fluoride improves MS polymer film's acid resistance. It protects the coated area from demineralization by dietary acids.

Inhibits acid erosion

● Tooth surface after
Gel Desensitizer application




● Protected area after soaked
in a soda for 10 min.



Protection from demineralization by lactic acid

(Application for 30 sec.)





The amount of calcium ion released by the acid challenge (lactate pH 4.5, for 2 hours) from the dentin pretreated with 3 different desensitizers.


Gel Desensitizer inhibited the demineralization as well as MS Coat F. Adding sodium fluoride contributed to the improved acid resistance of MS polymer film.

Decreased hypersensitivity at tooth whitening

Gel Desensitizer can reduce the hypersensitivity after tooth whitening.



Application of Gel Desensitizer before the whitening will not alter the whitening results.



The gel will not affect the future whitening treatments. The outcome of desensitized enamel will be as bright as the non-treated enamel.