

Super-Bond C&B

THE MOST REFERENCED GLUE
MILLIONS OF CLINICAL USES SINCE 1983
A TIME -PROVEN CHEMISTRY

*The Plus of
Super-Bond C&B*

*28 years of clinical success
all over the world !!*

Well known in main
universities



The only one

- with visco elasticity
- needs humidity for best curing
- almost complete curing (99 %)

THE KEY INDICATIONS OF SUPERBOND

BONDED BRIDGES

key points

- outstanding bonding values on alloys
- unrivaled technic



replacement of lateral front teeth



the natural bonding bridge :
the extracted crown is directly
bonded to proximal teeth

SPLINTING MOBILE TEETH

key points

- the splint is flexible !!
- very easy to use ,very esthetic result
- no secondary caries
- great after orthodontic treatment



INLAYS CORE CEMENTATION

key points

- self curing process in the root canal
- protection of the root due to absorption of constraints



INTRA ORAL REPAIR FOR CERAMICS

key points

- Porcelain Liner M + Superbond can repair any ceramic with any composite
- 15 years of experience with this association of 2 products



A VERY SIMPLE PRODUCT TO USE !

ONLY 1 mix

No bonding system necessary

A complete cure is achieved even in case of moisture or blood

Super-Bond Porcelain Liner M

Special primer
for ceramics



Porcelain Liner M is a coupling agent containing the famous 4 Meta monomer (from Superbond C&B) which should be used as a Primer in combination with Super Bond C&B.

Super Bond C&B
+ Porcelain Liner M :
20 years of ceramic bonding !!!

Clinical indications

- intra oral repair of fractured ceramics (PFM, All ceramics)
- bonding of inlays, onlays, veneers
- bonding of brackets on ceramic crown (orthodontics)

SEP isolating primer

A new blue primer to be applied on all external surfaces.

SEP will make the excess removal of any bonding so quick and easy Absolutely perfect to protect extrados and prevent proximal fillings from being touched by the bonding system.

● SEP is not fat

● Sep is hydrosoluble : you rince it in a second

Presentation :
5 ml bottle



Apply where you want : its blue and easy to check

V-PRIMER

The special Primer for
precious/semi precious alloys



Based on Dithiol Triazin, V Primer should be applied on sand-blasted Gold surfaces.

Apply, Dry and use Super Bond C&B .

Presentation :
3 ml bottle